


SEP15 Rabbit mAb

ABZ1851

Product Overview

Name	SEP15 Rabbit mAb
Accession(Primary)	Q9ERR7
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 18 kDa; Observed MW: 18 kDa
Recommended Dilution	WB-1:1000-1:5000

Target information (Q9ERR7)

Gene Name Selenof**About Gene Symbol**

Western blot analysis of 44089 in HeLa lysates using SEP15 antibody.

