

SMC1A Mouse mAb

ABZ4273


Product Overview

Name	SMC1A Mouse mAb
Accession(Primary)	Q14683
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human
Calculated MW	Calculated MW: 143 kDa; Observed MW: 143 kDa
Recommended Dilution	WB-1:500-1:1000; IF-1:50-1:200

Target information (Q14683)

Gene Name SMC1A

About Gene Symbol



Western blot analysis of SMC1A(Cterm) in K562, Jurkat, C6, 3T3 and HeLa lysates using SMC1A (Nterminus) antibody.

