


SPON2 Rabbit mAb

ABZ1870

Product Overview

Name	SPON2 Rabbit mAb
Accession(Primary)	Q9BUD6
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 36 kDa; Observed MW: 36 kDa
Recommended Dilution	WB-1:1000-1:5000

Target information (Q9BUD6)

Gene Name SPON2**About Gene Symbol**

Western blot analysis of SPON2 in K562, C6, 3T3, HeLa lysates using SPON2 antibody.

