

## STAT1 (9D6) Mouse mAb


ABZ3969

### Product Overview


|                             |                       |
|-----------------------------|-----------------------|
| <b>Name</b>                 | STAT1 (9D6) Mouse mAb |
| <b>Accession(Primary)</b>   | P42224                |
| <b>Host Species</b>         | Mouse                 |
| <b>Clonality</b>            | Monoclonal            |
| <b>Isotype</b>              | IgG1                  |
| <b>Species Reactivity</b>   | Human, Mouse, Rat     |
| <b>Recommended Dilution</b> | IHC-1:50-1:100        |

### Target information (P42224)

|                          |                                                                                                                                                                                                                                                                    |
|--------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Gene Name</b>         | STAT1                                                                                                                                                                                                                                                              |
| <b>About Gene Symbol</b> | STAT1 (Signal Transducer And Activator Of Transcription 1)<br>is a cytokine and immune signaling protein involved in immune signaling. It is an important research target in immunology, inflammation, and cancer research and is widely studied to elucidate cell |



Immunohistochemistry analysis of paraffin-embedded Human Breast Carcinoma Tissue using STAT1 (9D6) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunohistochemical analysis of paraffin-embedded Human tonsils using STAT1 (9D6) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.