

## STRA8 Rabbit pAb

ABZ1626


### Product Overview

|                             |                                            |
|-----------------------------|--------------------------------------------|
| <b>Name</b>                 | STRA8 Rabbit pAb                           |
| <b>Accession(Primary)</b>   | Q7Z7C7                                     |
| <b>Host Species</b>         | Rabbit                                     |
| <b>Clonality</b>            | Polyclonal                                 |
| <b>Isotype</b>              | IgG                                        |
| <b>Species Reactivity</b>   | Human, Mouse                               |
| <b>Calculated MW</b>        | Calculated MW: 37 kDa; Observed MW: 66 kDa |
| <b>Recommended Dilution</b> | WB-1:500-1:2000; IHC-1:50-1:300            |

### Target information (Q7Z7C7)

**Gene Name** STRA8

**About Gene Symbol**



Western blot analysis of STRA8 in NIH3T3, 293 lysates using STRA8 antibody.

