

STUB1 Rabbit mAb


ABZ1752

Product Overview


Name	STUB1 Rabbit mAb
Accession(Primary)	Q9UNE7
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 35 kDa; Observed MW: 35 kDa
Recommended Dilution	WB-1:1000-1:5000; ICC/IF-1:50-1:200; IP-1:20-1:50

Target information (Q9UNE7)

Gene Name	STUB1
About Gene Symbol	STUB1 (STIP1 Homology And U-Box Containing Protein 1) is a ubiquitin E3 ligase involved in ubiquitin-proteasome pathway and adaptive immunity. It is an important research target in cancer research and is widely studied to elucidate cellular mechanisms and



Western blot analysis of STUB1 in K562, rat Brain, C6, 3T3, HeLa lysates using STUB1 antibody.



Immunocytochemistry analysis of STUB1 (green) in HeLa using STUB1 antibody, and DAPI (blue)