

# Sodium/Hydrogen Exchanger 1/NHE-1 Rabbit mAb

ABZ5189

## Product Overview

|                             |                                              |
|-----------------------------|----------------------------------------------|
| <b>Name</b>                 | Sodium/Hydrogen Exchanger 1/NHE-1 Rabbit mAb |
| <b>Accession(Primary)</b>   | P19634                                       |
| <b>Host Species</b>         | Rabbit                                       |
| <b>Clonality</b>            | Monoclonal                                   |
| <b>Isotype</b>              | IgG                                          |
| <b>Species Reactivity</b>   | Human                                        |
| <b>Calculated MW</b>        | Calculated MW: 91 kDa; Observed MW: 130 kDa  |
| <b>Recommended Dilution</b> | WB-1:1000-1:5000                             |

## Target information (P19634)

**Gene Name** SLC9A1

**About Gene Symbol**



