

# TATA binding protein TBP Mouse mAb

ABZ4084


## Product Overview

|                             |                                    |
|-----------------------------|------------------------------------|
| <b>Name</b>                 | TATA binding protein TBP Mouse mAb |
| <b>Accession(Primary)</b>   | P20226                             |
| <b>Host Species</b>         | Mouse                              |
| <b>Clonality</b>            | Monoclonal                         |
| <b>Isotype</b>              | IgG1                               |
| <b>Species Reactivity</b>   | Human, Mouse, Rat                  |
| <b>Recommended Dilution</b> | IHC-1:50-1:100                     |

## Target information (P20226)

**Gene Name** TBP

### About Gene Symbol



Immunohistochemistry analysis of paraffin-embedded Human Breast Carcinoma using TATA Box Binding Protein (10C1) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.