

## TRIB3 Rabbit mAb


ABZ1896

### Product Overview

|                             |                                            |
|-----------------------------|--------------------------------------------|
| <b>Name</b>                 | TRIB3 Rabbit mAb                           |
| <b>Accession(Primary)</b>   | Q96RU7                                     |
| <b>Host Species</b>         | Rabbit                                     |
| <b>Clonality</b>            | Monoclonal                                 |
| <b>Isotype</b>              | IgG                                        |
| <b>Species Reactivity</b>   | Human                                      |
| <b>Calculated MW</b>        | Calculated MW: 40 kDa; Observed MW: 40 kDa |
| <b>Recommended Dilution</b> | WB-1:1000-1:5000; IP-1:10                  |

### Target information (Q96RU7)

|                          |                                                                                                                                                                                                                                                                                                                  |
|--------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Gene Name</b>         | TRIB3                                                                                                                                                                                                                                                                                                            |
| <b>About Gene Symbol</b> | TRIB3 (Tribbles Pseudokinase 3) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets. |



Western blot analysis of TRIB3 in K562, HeLa  
lysates using TRIB3 antibody.