

## USP10 Rabbit mAb


ABZ2122

### Product Overview


|                             |                                                                           |
|-----------------------------|---------------------------------------------------------------------------|
| <b>Name</b>                 | USP10 Rabbit mAb                                                          |
| <b>Accession(Primary)</b>   | Q14694                                                                    |
| <b>Host Species</b>         | Rabbit                                                                    |
| <b>Clonality</b>            | Monoclonal                                                                |
| <b>Isotype</b>              | IgG                                                                       |
| <b>Species Reactivity</b>   | Human                                                                     |
| <b>Calculated MW</b>        | Calculated MW: 87 kDa; Observed MW: 110 kDa                               |
| <b>Recommended Dilution</b> | WB-1:1000-1:5000; IHC-P-1:200-1:500;<br>ICC/IF-1:200-1:500; IP-1:10-1:100 |

### Target information (Q14694)


|                          |       |
|--------------------------|-------|
| <b>Gene Name</b>         | USP10 |
| <b>About Gene Symbol</b> |       |



Western blot analysis of USP10 in HeLa lysates using USP10 antibody.



Immunocytochemistry analysis of USP10(green) in HeLa using USP10 antibody, and DAPI(blue)



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using USP10 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.