

USP7 Rabbit mAb


ABZ1935

Product Overview


Name	USP7 Rabbit mAb
Accession(Primary)	Q93009
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 128 kDa; Observed MW: 128 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC:1:50-1:100; ICC/IF-1:100-1:200; FC-1:20-1:50

Target information (Q93009)

Gene Name	USP7
About Gene Symbol	USP7 (Ubiquitin Specific Peptidase 7) is a deubiquitinase involved in ubiquitin-proteasome pathway. It is an important research target in infectious disease and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of USP7 in C6 lysates using USP7 antibody.



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using HAUSP/USP7 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.