

Ubiquitin Rabbit mAb



ABZ3218


Product Overview

Name	Ubiquitin Rabbit mAb
Accession(Primary)	P0CG47
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 26 kDa; Observed MW: 8 kDa
Recommended Dilution	WB-1:1000-1:20000; IHC-P-1:100-1:200; ICC/IF-1:200-1:500; FC-1:100-1:1000

Target information (P0CG47)

Gene Name	UBB
About Gene Symbol	





Western blot analysis of Ubiquitin in 293 lysates
using Ubiquitin antibody