

CMV Pp65

Cytomegalo Virus Pp65 (UL83) Recombinant

VAG0058

Product Overview

Name

CMV Pp65

Description

Cytomegalo Virus Pp65 (UL83) Recombinant

References for protein:

Title: Cytomegalovirus-Infected Human Endothelial Cells Can Stimulate Allogeneic CD4+ Memory T Cells by Releasing Antigenic Exosomes Publication: Endocrinology 150.7 (2009): 3252-3258. Link: CMV Pp65 prospec publication