

CMV Pp65, 561 a.a.

Cytomegalo Virus Pp65(UL83), 561 a.a.Recombinant

VAG0059

Product Overview

Name	CMV Pp65, 561 a.a.
Immunogen	Anti-human C-Myc mAb, is derived from hybridization of mouse F0 myeloma cells with spleen cells from BALB/c mice immunized with a recombinant human C-Myc protein 410-419 amino acids purified from E. coli.
Description	Cytomegalo Virus Pp65(UL83), 561 a.a.Recombinant

References for protein:

Title: Cytomegalovirus-Infected Human Endothelial Cells Can Stimulate Allogeneic CD4+ Memory T Cells by Releasing Antigenic Exosomes Publication: Endocrinology 150.7 (2009): 3252-3258. Link: CMV Pp65.prospec publication