

HAVCR2 Mouse

Hepatitis A Virus Cellular Receptor 2 Mouse Recombinant

VAG0185

Product Overview

Name	HAVCR2 Mouse
Accession(Primary)	Q8VIM0
Immunogen	Anti-human HAX1 mAb is derived from hybridization of mouse F0 myeloma cells with spleen cells from BALB/c mice immunized with recombinant human HAX1 amino acids 1-279 purified from E. coli.
Description	Hepatitis A Virus Cellular Receptor 2 Mouse Recombinant

Target information (Q8VIM0)

Gene Name

About Gene Symbol