

HIV-1 p24 gag

HIV-1 p24 gag Recombinant

VAG0326

Product Overview

Name	HIV-1 p24 gag
Description	HIV-1 p24 gag Recombinant

References for protein:

Title: Polyfunctional CD4+ and CD8+ T Cell Responses to Tuberculosis Antigens in HIV-1-Infected Patients before and after Anti-Retroviral Treatment Publication: The Journal of Immunology 184.11 (2010): 6537-6544. Link: [HIV-1 p24 gag prospec publication](#)