

SARS Spike (408-470, 540-573)

SARS-Associated Coronavirus Spike (408-470, 540-573 a.a.), Recombinant

VAG0456

Product Overview

Name	SARS Spike (408-470, 540-573)
Accession(Primary)	P0C6X7.1
Immunogen	The antibody was developed using a synthetic peptide from the Spike S2 glycoprotein for the Human SARS coronavirus (Genbank accession number NP_828851.1) corresponding to amino acids 1124-1140.
Description	SARS-Associated Coronavirus Spike (408-470, 540-573 a.a.), Recombinant

Target information (P0C6X7.1)

Gene Name

About Gene Symbol