

HIV-1 gp120 CM

HIV-1 gp120 CM Envelope Recombinant

VAG0535

Product Overview

Name	HIV-1 gp120 CM
Accession(Primary)	L03703
Description	HIV-1 gp120 CM Envelope Recombinant

Target information (L03703)

Gene Name

About Gene Symbol

References for protein:

Title :Slit2/Robo4 Signaling Modulates HIV-1 gp120-Induced Lymphatic Hyperpermeability. Publication :Zhang X, Yu J, Kuzontkoski PM, Zhu W, Li DY, et al. (2012) Slit2/Robo4 Signaling Modulates HIV-1 gp120-Induced Lymphatic Hyperpermeability. PLoS Pathog 8(1): e1002461. doi:10.1371/journal.ppat.1002461 Link : <http://www.plospathogens.org/article/info%3Adoi%2F10.1371%2Fjournal.ppat.1002461>