

New York 3571/2009

Hemagglutinin-Influenza A Virus H1N1 New York 3571/2009 Recombinant

VAG0543

Product Overview

Name	New York 3571/2009
Immunogen	Anti-human NKp44 mAb, is derived from hybridization of mouse SP2/0 myeloma cells with spleen cells from BALB/c mice immunized with recombinant human NKp44 amino acids 19-130 purified from E. coli.
Description	Hemagglutinin-Influenza A Virus H1N1 New York 3571/2009 Recombinant