


XRCC1 Mouse mAb

ABZ4191

Product Overview

Name	XRCC1 Mouse mAb
Accession(Primary)	P18887
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human
Calculated MW	Calculated MW: 69 kDa; Observed MW: 82 kDa
Recommended Dilution	WB-1:500-1:1000

Target information (P18887)

Gene Name XRCC1**About Gene Symbol**

Western blot analysis of XRCC1 in Raji, K562, HeLa, Jurkat and THP-1 lysates using XRCC1 antibody.

