


ZFP106 Rabbit pAb

ABZ1474


Product Overview

Name	ZFP106 Rabbit pAb
Accession(Primary)	Q9H2Y7
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Cat
Calculated MW	Calculated MW: 209 kDa; Observed MW: 209 kDa
Recommended Dilution	WB-1:500-1:2000; ELISA-1:10000

Target information (Q9H2Y7)

Gene Name ZNF106**About Gene Symbol**

Western blot analysis of ZFP106 in COS7 lysates using ZFP106 antibody.



Western blot analysis of ZFP106 in COS7 lysates
using ZFP106 antibody.