

# alpha 1 Fetoprotein Rabbit pAb


ABZ1515

## Product Overview


<b>Name</b>	alpha 1 Fetoprotein Rabbit pAb
<b>Catalog #</b>	ABZ1515
<b>Clonality</b>	Polyclonal
<b>Accession(Primary)</b>	P02771
<b>Application Note (Approx.)</b>	WB1:500-1:2000 IHC1:100-1:300
<b>Precautions</b>	

## Target information(P02771)

**Synonyms****Gene ID** 174**Other Names****Function****Cellular location****Note**




Western blot analysis of alpha 1 Fetoprotein in mouse liver, mouse serum lysates using alpha 1 Fetoprotein antibody.



Immunohistochemistry analysis of paraffin-embedded rat brain using alpha 1 Fetoprotein antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunohistochemistry analysis of paraffin-embedded mouse brain using alpha 1 Fetoprotein antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunohistochemistry analysis of paraffin-embedded mouse brain using alpha 1 Fetoprotein antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.