

beta Arrestin 1 (6A1) Mouse mAb


ABZ3955

Product Overview

Name	beta Arrestin 1 (6A1) Mouse mAb
Accession(Primary)	P49407
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG2b
Species Reactivity	Human
Calculated MW	Calculated MW: 47 kDa; Observed MW: 51 kDa
Recommended Dilution	WB-1:500-1:1000

Target information (P49407)

Gene Name	ARRB1
About Gene Symbol	ARRB1 (Arrestin Beta 1) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of betaArrestin 1 in K562 and HeLa lysates using betaArrestin 1 antibody.