

eNOS (2E8) Mouse mAb

ABZ3989


Product Overview

Name	eNOS (2E8) Mouse mAb
Accession(Primary)	P29474
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human, Mouse, Rat, Rabbit
Calculated MW	Calculated MW: 133 kDa; Observed MW: 133 kDa
Recommended Dilution	WB-1:500-1:2000

Target information (P29474)

Gene Name NOS3

About Gene Symbol



Western blot analysis of eNOS (2E8) in rat Heart, HUVEC, Jurkat lysates using eNOS (2E8) antibody.

